


## Background

- The **Mayfield Head Clamp** (Figure 1) is a medical device that utilizes **3 pins** (Figure 2) to **hold the head rigid** during neurosurgery.
- A **pressure pin** is currently used to **communicate the pressure** exerted on the skull by the prongs.



## Problem

The **current method** of communicating pressure is **inaccurate, inefficient, and difficult to read**. This can lead to **complications** during surgery, including skull fractures, impaling of the brain, lacerations from the head slipping and superficial abrasions.

## Objectives

To **increase patient safety** during neurosurgery, the device should be designed to quickly **report the PSI exerted on the skull by the Mayfield Head Clamp** to neurosurgeons.

The device must be:

- Sterilizable** and compatible with hospital cleaning standards to protect patient safety
- Integrable** with the current clamp and not obstruct standard practices
- Capable of **withstanding** the safe **pressure threshold** applied during neurosurgery

## Results




Figure 3: Final prototype



Figure 4:  
Aluminum pins  
with elongated divot




Figure 5:  
Calibration Curve

**Final Prototype**  
(Figure 3) **Includes:**

- Aluminum** pins with **elongated divot** to accommodate a strain gauge (Figure 4)
- Epoxy-attached **strain gauge** to measure deformation
- Calibration curve** to convert from deformation to PSI readings (Figure 5)
- Screen & LED** to display pressure (in PSI) (Figure 6)




Figure 6:  
Screen/LED integration

## Design Criteria and Testing

Design Criteria	Target Value	Test	Results	
Accuracy	Pressure readings <b>within ±5 PSI, 95%</b> of the time	Comparing reported PSI to PSI indicated by torque screw on clamp	Passed <b>5/5</b> calibration curve tests	✓
Intuitive to Use	Pressure readings take <b>&lt;1 minute</b>	User given pre-tightened clamp and asked to return PSI	Readings took <b>~10 seconds</b>	✓
Pressure Threshold	Strain gauge can <b>withstand up to 120 PSI</b> , 80 PSI minimum	Recording voltage readings at 60, 80, 100 and 120 PSI	Dynamic range of strain gauge is <b>150 PSI</b>	✓
Durability	<b>5 year</b> lifespan	Research and lifespan calculation	Pins are <b>replaceable</b> and electronics last <b>10+ years</b>	✓
Sterilizability	<b>Withstand 15 sterilization cycles</b> at 134°C	N/A due to available resources	TBD	

## Conclusion

Our team modified the Mayfield Head Clamp pins to allow for the ability to **efficiently** and **accurately** **inform** neurosurgeons of the **pressure** exerted **on the skull** by the clamp during surgery.

Future enhancements/plans may include:

- Increasing prototype fidelity by moving from aluminum to stainless steel
- Further calibration and accuracy testing
- Sterilizability testing

## Acknowledgements

Our team would like to thank the following for their continued support throughout this project.

**Dr. Sonia Bansal**, Professor  
**Dr. Micheal Rizk**, Professor  
**Dr. Eli Johnson**, Client  
**Kevin Shores**, Technical Mentor  
**Alyssa Ramirez**, Teaching Assistant  
**Erin O'Rourke**, Writing Consultant



## Background

- The **Mayfield Head Clamp** (MHC) (Figure 1) is a medical device that utilizes **3 pins** (Figure 2) to **hold the head rigid** during neurosurgery.
- An **inefficient pressure pin** is currently used to **communicate the pressure** exerted on the skull by the prongs.



## Problem

The **current method** of communicating pressure is **inaccurate, inefficient, and difficult to read**. This can lead to **complications** during surgery, including skull fractures, impaling of the brain, lacerations from the head slipping and superficial abrasions.

## Objective

To **increase patient safety** during neurosurgery, the device should be designed to quickly **report the PSI** exerted on the skull by the Mayfield Head Clamp to neurosurgeons.

## Design Criteria and Testing Plan

Design Criteria	Target Value	Test
Accuracy / Sensitivity	Pressure readings <b>within ±5 PSI</b> , <b>95%</b> of the time	Comparing reported PSI to PSI indicated by torque screw on clamp. Additional verification using material properties.
Integratable	Pins can fit into a <b>variety of devices</b>	Testing usability across multiple brands of the Mayfield Head Clamp.
Intuitive to Use	<b>Digital</b> pressure readings take <b>&lt;1 minute</b>	User given pre-tightened clamp and asked to return PSI
Pressure Threshold	Strain gauge can <b>withstand up to 120 PSI</b> , 80 PSI minimum	Recording voltage readings at 60, 80, 100 and 120 PSI
Sterilizability	<b>Single use</b> pins	N/A

## Results





Figure 5:  
Screen/LED display integration




Figure 6:  
60mm x 60mm PCB to digital display integration

## Feature Comparison

Feature	Iteration 1	Iteration 2
Aluminum pin	✓	✓
Elongated divot	✓	✓
Flattened divot		✓
Epoxy-attached Strain Gauge	✓	✓
Flush fit		✓
Screen and LED display	✓	
60x60 mm PCB to digital display		✓

### Pin Modification:

- Flattened divot** increases **durability** of the pin to strain gauge connection
- Flush fit** allows for **increased accuracy** of the strain gauge
- More precise measurements allow for a **smoother fit** into the clamp

### Display Modification:

- Personalized PCB** allows for easier modification to accommodate different strain gauges
- PCB can connect to various devices**
- Decreased size** reduces visual disorder and makes display **more intuitive**

## Conclusion

Our team modified the Mayfield Head Clamp pins to allow for the ability to **efficiently** and **accurately inform** neurosurgeons of the **pressure** exerted **on the skull** by the clamp during surgery.

Future enhancements/plans may include:

- Increasing prototype fidelity by moving from aluminum to stainless steel
- Further calibration and accuracy testing
- Sterilizability testing

## Acknowledgements

Our team would like to thank the following for their continued support throughout this project.

Dr. Ibrahim Mohedas , Professor  
Dr. Ann Saterbak, Professor  
Dr. Eli Johnson, Client