



Pressure Sensor for Mayfield Head Clamp

Saagar Arya, Anna Brusoe, Lindsay Chetkof, Amanda Li Duke University, Duke Neurosurgery



Background

- The Mayfield Head
 Clamp (Figure 1) is a
 medical device that utilizes
 3 pins (Figure 2) to hold
 the head rigid during
 neurosurgery.
- A pressure pin is currently used to communicate the pressure exerted on the skull by the prongs.



Mayfield Head

Clamp pin

Problem

The current method of communicating pressure is inaccurate, inefficient, and difficult to read. This can lead to complications during surgery, including skull fractures, impaling of the brain, lacerations from the head slipping and superficial abrasions.

Objectives

To increase patient safety during neurosurgery, the device should be designed to quickly report the PSI exerted on the skull by the Mayfield Head Clamp to neurosurgeons.

The device must be:

- Sterilizable and compatible with hospital cleaning standards to protect patient safety
- Integrable with the current clamp and not obstruct standard practices
- Capable of withstanding the safe pressure threshold applied during neurosurgery

Results



Figure 3: Final prototype




Figure 5:Calibration Curve

Final Prototype (Figure 3) Includes:

- Aluminum pins with elongated divot to accommodate a strain gauge (Figure 4)
- Epoxy-attached strain gauge to measure deformation
- Calibration curve to convert from deformation to PSI readings (Figure 5)
- Screen & LED to display pressure (in PSI) (Figure 6)




Figure 6:
Screen/LED integration

Design Criteria and Testing

Design Criteria	Target Value	Test	Results	
Accuracy	Pressure readings within ±5 PSI, 95% of the time	Comparing reported PSI to PSI indicated by torque screw on clamp	Passed 5/5 calibration curve tests	
Intuitive to Use	Pressure readings take <1 minute	User given pre-tightened clamp and asked to return PSI	Readings took ~10 seconds	
Pressure Threshold	Strain gauge can withstand up to 120 PSI, 80 PSI minimum	Recording voltage readings at 60, 80, 100 and 120 PSI	Dynamic range of strain gauge is 150 PSI	
Durability	5 year lifespan	Research and lifespan calculation	Pins are replaceable and electronics last 10+ years	
Sterilizability	Withstand 15 sterilization cycles at 134°C	N/A due to available resources	TBD	

Aluminum pins

with elongated divot

Conclusion

Our team modified the Mayfield Head Clamp pins to allow for the ability to **efficiently** and **accurately inform** neurosurgeons of the **pressure** exerted **on the skull** by the clamp during surgery.

Future enhancements/plans may include:

- Increasing prototype fidelity by moving from aluminum to stainless steel
- Further calibration and accuracy testing
- Sterilizability testing

Acknowledgements

Our team would like to thank the following for their continued support throughout this project.

Dr. Sonia Bansal, Professor
Dr. Micheal Rizk, Professor
Dr. Eli Johnson, Client
Kevin Shores, Technical Mentor
Alyssa Ramirez, Teaching Assistant
Erin O'Rourke, Writing Consultant